



Shimadzu Diagnostics Corporation

Culture Medium for Regenerative and Cell Therapy Research

Liquid Medium

02	MSF	06	NS-A2, DCO-K
03	MSH	07	PSE8
04	MSH-EV	08	Contract Media Manufacturing Service
05			Culture Medium Component Analysis Service

Features

- Medium for the culture of human mesenchymal stromal cells**

When MSF-BM and Supplement A are used in combination, human mesenchymal stromal cells (MSC) can be cultured at a high proliferation rate while maintaining their undifferentiated state, compared to classical media with fetal bovine serum (FBS)

- Choice of culture additives**

Supplement A is a low serum (2% final concentration) supplement developed for MSC expansion culture (FBS is γ -ray treated) Alternatively, combination with the human platelet-derived cell culture additive UltraGRO can increase proliferation


- Certified as eligible for regenerative medical products**

MSF-BM has obtained a certificate of material eligibility for regenerative medical products from PMDA (Pharmaceuticals and Medical Devices Agency)

Adipose-derived MSC proliferation study using MSF-Supplement A




MSF-BM + Supplement A provides better support for proliferation than α MEM with 20% FBS.



Adipose-derived MSC proliferation study using UltraGRO

5% UltraGRO was added to each medium, and the proliferation rates of adipose-derived MSC compared.

Higher proliferation of MSC was obtained with MSF-BM than with the classical media, D-MEM and α MEM.



Features

- **Xeno-free medium**

It does not contain fetal bovine serum (FBS) or human platelet lysate

The protein ingredients consist solely of recombinant proteins expressed in microorganisms, and pharmaceutical grade human serum-derived albumin

- **Good support for proliferation**

Higher proliferation compared to α MEM with 20% FBS

- **Maintains an undifferentiated state**

Proliferates human mesenchymal stromal cells (MSC) while maintaining their differentiation potential



Marker expression and differentiation potential into adipocytes, osteoblasts, and chondrocytes confirmed



MSC proliferation tests and surface markers




Higher proliferation of human adipose tissue-derived MSC was obtained in MSH medium than in α MEM with 20% FBS.

In addition, positivity for MSC markers CD73, CD90, and CD105 was confirmed.



Confirmation of MSC differentiation potential

Human adipose tissue-derived MSC cultured in MSH medium for six passages maintained the ability to differentiate into adipocytes, osteoblasts, and chondrocytes.



Features

- Chemically defined medium**

Can be used as a medium for the proliferation of human mesenchymal stromal cells (MSC) and exosome production

- Improved exosome purification**

MSH-EV is low in particulate contamination*, enabling the use of ultracentrifugation and ultrafiltration
Reduces contamination by particles derived from the culture medium during exosome purification

- Exosomes can be recovered from the supernatant during expansion culture**



A series of processes — from MSC proliferation to exosome purification — can be performed in a single medium

* Lower than the detection limit by Nanoparticle Tracking Analysis (model: NanoSight)






MSC proliferation potential

MSH-EV medium could support the proliferation of human adipose tissue-derived MSC (hAd-MSC) at an equal or higher level than α MEM with 20% FBS.




Confirmation of MSC differentiation potential

When hAd-MSC was cultured in MSH-EV medium for approximately 1 month, it maintained the ability to differentiate into adipocytes, osteoblasts, and chondrocytes.




MSC exosome production capacity (as medium for expansion culture)

Culture supernatants from hAd-MSC expanded and cultured in MSH-EV medium yielded as many exosomes as those using α MEM + 20% FBS (exosomes positive for phosphatidylserine and CD63 by ELISA).



Exosome recovery and particle size distribution

MSH-EV medium yielded more exosomes than α MEM or competitor's exosome recovery media (exosomes positive for Phosphatidylserine and CD63 by ELISA).



List of Related Products

Product Name	Code	Capacity	Storage temperature	Purpose
MSF-BM	66023	500 mL	2–10°C	Human mesenchymal stromal cell culture medium
MSF-Supplement A	66024	10 mL	below –20°C	MSF-BM supplement (Hypoglycemia)
MSH-BM	66213	500 mL	2–10°C	Serum-free, human mesenchymal stromal cell culture medium
MSH-Supplement A	66214	10 mL (for 500 mL)	below –20°C	MSH-BM supplement (Xenofree)
MSH-EV-Supplement A	66216	0.2 mL (for 100 mL)	below –20°C	MSH-BM supplement (for EV production and recovery)
DCO-K	66017	1,000 mL	2–10°C	Serum-free culture medium for immunotherapy research
NS-A2	66001	1,000 mL	2–10°C	Serum-free culture medium for lymphocytes

This product line consists of research reagents.

Serum-free culture medium for lymphocytes

NS-A2

Research Reagent

NS-A2 is a serum-free culture medium developed for immunotherapy applications.

- Supports high proliferation of human peripheral blood T cells.
- It does not use animal-derived protein ingredients.
- CD3 antibodies and other agents support the high proliferation of activated peripheral blood T lymphocytes.
- No protein ingredients other than pharmaceutical-grade human serum albumin and recombinant human insulin are used.
- Testing for sterility, endotoxins, and mycoplasma negativity in accordance with the Japanese Pharmacopoeia are conducted in addition to physicochemical tests.



Literature

Long-term eradication of extranodal natural killer/T-cell lymphoma, nasal type, by induced pluripotent stem cell-derived Epstein-Barr virus-specific rejuvenated T cells *in vivo*. Miki Ando., et al. *Haematologica* 2020 Volume 105(3):796-807.

<https://doi.org/10.3324/haematol.2019.223511>

Dual-antigen targeted iPSC-derived chimeric antigen receptor-T cell therapy for refractory lymphoma. Sakiko Harada., et al. *Molecular therapy* Vol. 30 No 2: 534-549.

<https://doi.org/10.1016/j.ymthe.2021.10.006>

iPSC-derived hypoimmunogenic tissue resident memory T cells mediate robust anti-tumor activity against cervical cancer. Yoshiki Furukawa., et al. *Cell Rep Med*. 2023 Dec 8:101327.

<https://doi.org/10.1016/j.xcrm.2023.101327>

DCO-K

Research Reagent

DCO-K is a serum-free culture medium developed for immunotherapy research.

- Enables efficient differentiation from human monocytes into dendritic cells.
- Does not use animal-derived protein ingredients.
- Addition of UltraGRO, etc. improves yield and survival rate after culture.
- Testing for sterility, endotoxins, and mycoplasma negativity in accordance with the Japanese Pharmacopoeia are conducted in addition to physicochemical tests.
- The product has obtained a certificate of eligibility for regenerative medical product materials from PMDA (Pharmaceuticals and Medical Devices Agency).



Literature

Interferon- α -Induced Dendritic Cells Generated with Human Platelet Lysate Exhibit Elevated Antigen Presenting Ability to Cytotoxic T Lymphocytes. Ippei Date., et al. *Vaccines* 2021,9(1), 10.

<https://dx.doi.org/10.3390/vaccines9010010>

Clinical Trial on the Safety and Tolerability of Personalized Cancer Vaccines Using Human Platelet Lysate-Induced Antigen-Presenting Cells. Terutsugu Koya., et al. *Cancers* 2023, 15, 3627.

<https://doi.org/10.3390/cancers15143627>

Features

- **Chemically defined medium with a composition similar to Essential 8 (E8)**
Original composition developed in-house (manufactured in Japan).
Smooth switching of culture media from E8.
- **Does not use animal-derived ingredients**
All protein ingredients are recombinants expressed in microorganisms or plants.
- **No need to exchange media on weekend**
By using ingredients with high thermal stability, we have improved the ease of culture operations.

Recommended use

- **Passaging method:** clump passaging using EDTA as cell detachment medium is recommended
- **Scaffolding material:** Vitronectin is recommended

Recommended schedule for medium exchange and passaging

Mon.	Tue.	Wed.	Thu.	Fri.	Sat.	Sun.
Passaging	Medium exchange	—	Passaging	Medium exchange (double dose)	—	—


Evaluation data: proliferation support capacity

Evaluation method

- **Human iPS cells:** strain 802-3G (Reprocell)
- **Cell detachment medium:** PBS containing 0.5 mM EDTA (clump passaging)
- **Culture vessel:** Vitronectin-coated 6-well plate

- **Evaluation medium:** (ROCK inhibitor added to each medium only when passaging)
 1. PSE8-BM/Supplement A
 2. Competitor's E8 series medium (Weekend-free)
- **Cell counting:** Vi-CELL XR (Beckman Coulter)
- See "Recommended use" for medium exchange and passaging schedules


Colony form



Colony morphology of iPS cells cultured in PSE8 was similar to that of competitor's E8-based medium (Weekend-free)


* We have a track record of passaging for over 90 days.

Proliferation of iPS cells




Without an exchange of medium on the weekend, PSE8 could support proliferation at almost the same level as the competitor's E8-based medium (Weekend-free)

Evaluation data: maintenance of undifferentiation



Evaluation data: differentiation potential

Human iPS cells cultured in PSE8 for 21 passages maintained their differentiation potential



Evaluation method: Human iPS cells cultured in PSE8 for 21 passages were differentiated and stained with markers as follows

- **Ectoderm differentiation:** cultured in DMEM/F12+5% KSR for 2 weeks to form embryoid bodies (EB). The resulting EB were dispersion treated with EDTA, then adhered to gelatin-coated plates and incubated for an additional 5 days
- **Mesoderm differentiation:** Cultured in STEMdiff Trilineage Mesoderm Medium for 5 days
- **Endoderm differentiation:** Cultured in STEMdiff Trilineage Endoderm Medium for 5 days

Under development

NK-Cell, HEK293

Contract Media Manufacturing Service

We can provide cell culture media tailored to your specific requirements (Custom formulations, modified compositions, packaging/volume, etc.). We offer meticulous support through prior consultation.

Features

- Operated under ISO 9001 Quality Management System
- Facility operation based on PIC/S GMP Annex 1
- Only raw materials compliant with STANDARDS FOR BIOLOGICAL RAW MATERIALS
- No use of highly bioactive substances such as antibiotics
- Production is possible in accordance with raw material standards for regenerative medicine products (powder and liquid culture medium)

Regenerative Medicine Culture Medium Plant (GMP-Compliant Manufacturing)

Powdered Medium : Up to 10 kg/batch

Liquid Medium : Up to 500 L/batch

Research Reagent Grade Manufacturing Available

Paid Samples for Evaluation Also Available

Inquire regarding manufacturing conditions and scale



Culture Medium Component Analysis Service

Analyzing 95 components in culture media using LCMS

Profiling 95 components (Medium components and secreted metabolites)

Amino Acid and Derivatives		Vitamins	Nucleic acid associated compounds	Others	Internal Standard
2-Aminoadipic acid	Histidine	4-Aminobenzoic acid	Adenine	2-Aminoethanol	2-Isopropylmalic acid
4-Aminobutyric acid	Isoleucine	Ascorbic acid	Adenosine	2-Ketoisovaleric acid	
4-Hydroxyproline	Kynurenine	Ascorbic acid 2-phosphate	Adenosine monophosphate	3-Methyl-2-oxovaleric acid	Sugars
5-Glutamylcysteine	Leucine	Biotin	Cytidine	4-Hydroxyphenyllactic acid	
5-Oxoproline	Lysine	Choline	Cytidine monophosphate	Citric acid	
Alanine	Methionine	Cyanocobalamin	Deoxycytidine	Ethylenediamine	
Alanyl-glutamine	Methionine sulfoxide	Ergocalciferol	Guanine	Fumaric acid	
Arginine	N-Acetylaspartic acid	Folic acid	Guanosine	Glyceric acid	
Asparagine	N-Acetylcysteine	Folinic acid	Guanosine monophosphate	Histamine	
Aspartic acid	Ornithine	Lipoic acid	Hypoxanthine	Isocitric acid	Gluconic acid
Citrulline	Oxidized glutathione	Niacinamide	Inosine	Lactic acid	Glucosamine
Cystathionine	Phenylalanine	Nicotinic acid	Thymidine	Malic acid	Hexose (Glucose)
Cysteine	Pipecolic acid	Pantothenic acid	Thymine	O-Phosphoethanolamine	Sucrose
Cystine	Proline	Pyridoxal	Uracil	Putrescine	Threonic acid
Glutamic acid	Serine	Pyridoxine	Uric acid	Pyruvic acid	
Glutamine	Threonine	Riboflavin	Uridine	Succinic acid	
Glutathione	Tryptophan	Tocopherol acetate	Xanthine	Penicillin G	
Glycine	Tyrosine		Xanthosine		
Glycyl-glutamine	Valine				

This product line consists of research reagents.

Copyright © Shimadzu Diagnostics Corporation. All rights reserved.

Shimadzu Diagnostics Corporation

E-mail: contact@sdcs.shimadzu.co.jp

TEL: +81 (0) 3 (5846) 5701

URL: <https://corp.sdc.shimadzu.co.jp/english/>



BO2603