

Case

CompactDry × Commercial Laboratory



Customers



Intertek Pakistan Pvt Ltd
accredited 17025 & ISO 9001 Karachi, Pakistan
<https://www.intertek.com/pakistan/>

Intertek Pakistan is a wholly owned subsidiary of Intertek Group PLC, U.K. They cater for commercial inspection and testing requirement of Food Division includes Microbiological testing, Shelf life, Residues, Analytical services and other chemical testing.

They are working with Multinational and national firms such as Nestle, Unilever, United Nation, PepsiCo, McDonald's, Coca Cola, National foods, Shan Foods, Asian foods, Young' s Food, Qarshi, EBM, Rahmat-e-Shireen, Metro, Haleeb Foods, Fauji Foods, Mitchell's, Shezan, Six B Foods, Bunny's & Punjab Food Authority.

Intertek Pakistan Lahore Lab/Office



Accreditation for Microbiological Laboratories

According to a guidance by Eurachem who is a network of organisations in Europe, having the objective of establishing a system for the international traceability of chemical measurements and the promotion of good quality practices,

In terms of Validation and Verification of test methods which is equivalent to ISO17025 7.2, it says that ""The following standards can assist laboratories in obtaining method validation data: ISO/TR 13843, **ISO16140** and ISO17994."" And, ""This validation data may be obtained through collaborative testing and from validation data submitted by the manufacturers and subjected to third party evaluation (e.g. AFNOR, **NordVal**, **Microval**, AOAC).""




Sources;
<https://www.eurachem.org/index.php>

ISO16140

ISO16140 is "Microbiology of the food chain- Method validation-." And, Part2 is "Protocol for the validation of alternative(proprietary) methods against a reference method."

INTERNATIONAL
STANDARD
ISO 16140-2

CompactDry



The validation and certification has been performed in accordance with EN ISO16140-2:2016 and the MicroVal Rules and Certification Scheme version 8.

NordVal International Certificate

Issued for: HyServe Compact Dry YM RAPID
NordVal No: 050
Approval date: 15 June 2018
Valid until: 15 June 2020

HyServe Compact Dry YM RAPID

Manufactured by: Nissui Pharmaceutical Co. Ltd, 3-23-9 Ueno, Taitoh-ku, Tokyo, 110-8736, Japan
Supplied by: HyServe GmbH & Co. KG, Hohenrainstr. 24, 52449 Lutzerath, Germany

fulfils the requirements of ISO 16140-2:2016 and the NordVal validation protocol 1. The reference method was ISO 21527-1:2009: Microbiology of food and animal feeding stuffs - Horizontal method for the enumeration of yeasts and moulds - colony count technique - Part 1: Products with a water activity greater than 0.05.

fulfils the requirement of ISO16140-2:2016 and the NordVal validation protocol 1.

Shimadzu Diagnostics Corporation

Scientific Laboratory Supplies

Address: Plot No. F14, Suit No. 301, Third Floor
S.I.T.E. Super Highway Phase-2 Karachi, Pakistan

Call Us: +92 300 2511883

**Local
Distributor**