


Case

CompactDry × nutritional milk powder



Customers



A leading Dairy manufacturer
in Vietnam

Milk powder products are, in general, more sensitive than liquid milk products, microbiologically. The manufacturing method of liquid milk products enables the aseptic or the sterilization, while milk powder products do not. Therefore, on the basis of the pasteurization, a strict microbiological control is required in all of the manufacturing process.


This customer adopts several international food safety system such as Haccp, ISO22000, FSSC22000, and provides both of domestic and overseas markets with the milk powder products in good quality.




Applications

This customer uses CompactDry implement the microbiology test in a laboratory. CompactDry is the ready to use culture media which is validated by the accrediting body such as AOAC, MicroVal and NordVal. You can not only save the time of media preparations, but also have more rapid detection time than the traditional culture method.


CompactDry "YM(Yeasts & Moulds) Rapid" is developed based on Potato dextrose agar media. The incubation time is from 3 days.




Total viable count



E.coli; blue Coliforms; red



YMR



Shimadzu Diagnostics Corporation

Local
Distributor



PACIFIC SCIENCE TECHNOLOGY EQUIPMENT CO.,LTD

Address: Floor 3rd, No. 382/17-19 Nguyen Thi Minh Khai str,
Ward 5, District 3, HCMC, VN

Web: <http://www.pacificlab.vn/> **E-Mail:** sales@pacificlab.vn